



Downloading HANDi Stroke Rx

Adam B. Schlichting

The process of downloading and installing HANDi Stroke Rx from the FERNE website has been designed to be as rapid and easy as possible. These directions will allow you to begin using HANDi Stroke Rx on your handheld computer in less than five minutes.

Downloading HANDi Stroke Rx

To download and install HANDi Stroke Rx, follow these simple instructions:

1. On your computer, open your web browsing program such as Netscape Navigator or Microsoft Internet Explorer.
2. Go to the FERNE website, <http://www.FERNE.org>
3. Click on the “Download the Stroke Management Program for Handhelds” button.
4. From the “Software” page, select the “NIHSS, Stroke Management, and t-PA Administration” link.
5. Select the type of computer system you are running (IBM/PC or Macintosh) from the Handheld Computer Stroke Program screen.
6. Input your name and email address on the “Handheld Computer Stroke Program Download” page, then click the “Submit” button.

Follow the operating system-specific instructions for the remainder of the installation process.

PC/IBM

7. RIGHT click the “Download” link, and select the “Save Target As...” option. Save the program to your computer desktop.
8. Double click on the “HANDi Stroke Rx” icon on your Desktop.
9. Follow the onscreen instructions.
10. HANDi Stroke Rx will be installed and ready to use next time you HotSync your Palm device!

Macintosh

7. Click the “Download” link and save the file to your Desktop.
8. Unzip HANDi Stroke Rx from your Desktop.
9. Open the Palm Desktop program on your computer.
10. Add each of the seven unzipped files to be installed using the Install Tool of the Palm Desktop software.
11. HANDi Stroke Rx will be installed and ready to use next time you HotSync your Palm device!

Questions regarding the program or its installation can be emailed to FERNE@FERNE.org